



MULTI- FUNCTIONAL VIAL OPENER

A Safe and Handy Vial Opening Assistant

This product is not classified as a medical device.



MULTI-FUNCTIONAL VIAL OPENER

Multi-functional Vial Opener:

Multi-functional vial opener was designed to be a safer tool to open various sizes of vials.

Not only providing full protection during operation, but also reducing the risk from stabbed by sharp stuff.


Feature:

- Attached blade for a quick removal of the crimped aluminum seals from the vial or bottle.
- Dual jaws for opening different sizes of cap.
- Safe & Easy to operate.
- Non-slip handle.
- Safety locking.



SPECIFICATIONS

Product	REF. No.	Size	Package	Description
Multi-functional Vial Opener	9001-9001	Universal	1 pc/case	Large jaw : For cap with O.D. 21mm-38mm. Small jaw : For cap with O.D. 18mm-26mm.



9001-9001



6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,
Danshuei Dist, New Taipei City. 251 Taiwan
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266
www.fortunemed.com Mail:info@fortunemed.com